



Digital Color Doppler Palm Ultrasound System

Model
SonoEye P6/SonoEye V6/SonoEye G6

USER MANUAL

Sep. 3, 2025

Rev. 1.0



CHISON Medical Technologies Co., Ltd.

We reserve the right to make changes to this manual without prior notice.

Notice d'information

Documentation utilisateur complète disponible sur demande

Vous consultez actuellement une version synthétique de la documentation relative à ce dispositif médical.

Seule une version allégée du manuel d'utilisation est mise à disposition sur notre site internet. Le manuel utilisateur complet peut être obtenu gratuitement sur simple demande.

Pour recevoir la version intégrale du manuel, nous vous invitons à compléter le formulaire de contact disponible sur notre site internet en précisant la référence du produit concerné.

Notre équipe vous transmettra la documentation complète dans les meilleurs délais, sous format électronique.

Pour toute question complémentaire, notre service client reste à votre disposition.

Important : Avant toute utilisation du dispositif, assurez-vous d'avoir pris connaissance des informations et consignes adaptées à votre usage.

La direction

Regulatory Requirement

CE 0197 This product conforms to the essential requirements of the Medical Device Regulation (EU) 2017/745. Accessories without the CE mark are not guaranteed to meet the Essential Requirements of the Medical Device Regulation.

This manual is a reference for the SonoEye P6/SonoEye V6/SonoEye G6. Please verify that you are using the latest revision of this document. If you need to know the latest revision, contact your distributor.

NOTE:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

TABLE OF CONTENTS

Chapter 1	Introduction	1
1.1	System Overview	1
1.2	Contact Information.....	1
Chapter 2	System Safety	2
2.1	Safety Overview.....	2
2.2	Electrical Safety	3
2.3	Labels.....	5
2.3.1	Warning Symbols.....	5
2.4	Patient Environmental Devices	7
2.5	Biological Safety	8
2.6	Scanning Patients and Education.....	9
2.7	Device instructions of cybersecurity controls.....	14
2.8	Insturctions to the intruded device	15
Chapter 3	Preparing the System for Use	17
3.1	Site Requirement	17
3.1.1	Operation Environmental Requirement.....	17
3.1.2	Electrical Requirements	17
3.1.3	Internal battery specifications	17
3.1.4	Hardware Requirements.....	18
3.2	Downloading and Installing the App	20
3.3	System Specifications	22
3.3.1	Console Overview	22
3.3.2	Physical Specifications.....	22
3.3.3	Image Modes	22
3.3.4	Type	23
3.3.5	System Configuration	24
3.4	System Positioning & Transporting.....	26
3.5	Powering the System	27
3.6	The System	27
Chapter 4	Control Panel	29

4.1	Overview of Display Area	29
4.2	Start a New Exam	29
4.3	Application	30
4.4	Switching Function	30
4.5	Switching Mode	31
4.6	Function Key Introduction	31
4.7	Image Parameter Area	35
4.8	Cine Control.....	35
Chapter 5 Imaging		36
5.1	Select Scan Mode	36
5.1.1	System Identify	36
5.1.2	Select Application	36
5.1.3	B mode	36
5.1.4	2B mode	36
5.1.5	FHI.....	36
5.1.6	B/M Mode.....	37
5.1.7	CFM Mode	37
5.1.8	CPA mode	37
5.1.9	DPD mode	37
5.1.10	PW mode	37
5.1.11	Biopsy Guide	38
5.1.12	SonoNeedle	38
5.1.13	SonoColor	38
5.1.14	SonoHelp	39
5.1.15	SonoCam	39
5.1.16	SonoCalc.....	39
5.1.17	SonoShare.....	40
5.1.18	SonoRemote.....	40
5.1.19	SonoBladder	41
5.1.20	SonoCalc.....	42
5.1.21	Full Screen.....	42
5.1.22	SonoZoom.....	42
5.2	Functional Description of Parameter Adjustment	42

5.3	Parameter Adjustment in PW mode.....	43
5.4	After Capturing the Image	44
5.4.1	Adding Annotation	44
5.4.2	Adding BodyMark.....	44
5.4.3	Adding Mark.....	45
Chapter 6 Measurement and Calculation.....		46
6.1	Fast measurement.....	46
6.1.1	Fast measurement in B Mode.....	46
6.1.2	Fast measurement in B/M Mode	47
6.1.3	Fast measurement in PW Mode	47
6.2	Measurement Package	48
6.2.1	Vessel Measurement in B mode	48
6.2.2	Abdomen Measurement in B mode.....	49
6.2.3	Cardiac Measurement in B mode.....	50
Chapter 7 Cine-Memory		52
7.1	Image Storage	52
7.2	Manual playback.....	52
7.3	Automatic playback.....	52
7.4	Cine Save	52
Chapter 8 Side Menu.....		53
8.1	Archive	53
8.2	EasyView	54
8.3	Report.....	55
8.4	Setting.....	56
8.4.1	General	56
8.4.2	Language	56
8.4.3	Button Configuration	56
8.4.4	Display Management.....	56
8.4.5	Reset.....	56
8.4.6	Preset	56
8.4.7	Optimize SonoEye Usage	57
8.4.8	Administrator	57
8.4.9	Firmware Update.....	57

8.4.10	Print Management.....	57
8.4.11	User Define Template Management.....	58
8.5	Save Preset	58
8.6	Export Destinations	58
8.7	Modality Worklist	58
8.8	SonoRemote	58
8.9	Tutorials	58
8.10	Demo.....	58
8.11	SonoShare.....	59
8.12	About	59
8.12.1	Function Management	59
8.12.2	Help and Feedback	59
8.12.3	Remote Service.....	59
8.12.4	Check for Upgrade	59
Chapter 9	System Maintenance.....	60
9.1	Cleaning.....	60
9.2	System Maintenance.....	60
9.3	Safety Check	61
9.4	Troubleshooting.....	61
9.5	Service Responsibility	62
Chapter 10	System.....	63
10.1	General Description.....	63
10.2	Care and Maintenance	63
10.2.1	Inspecting Systems.....	63
10.2.2	Preparation at the Point of Use	64
10.3	Cleaning and Disinfecting.....	64
10.4	System Operation Instructions.....	67
10.5	Service Responsibility	67
Appendix A:	THE INFORMATION OF EC REPRESENTATIVE	69
Appendix B:	ACOUSTIC OUTPUT REPORT TABLE	70
Appendix C:	TRANSDUCER MAXIMUM SURFACE TEMPERATURE.....	74
Appendix D:	MEASUREMENT RESULTS SUMMARY	75
Appendix E:	GUIDANCE AND MANUFACTURER'S DECLARATION.....	76

Appendix F: Pulsed Wave Doppler Velocity Measurement Results Summary..... 79

Chapter 1 Introduction

This manual contains necessary information for safe system operation.

Read and understand all instructions in this manual before operating the system. Always keep this manual with the equipment, and periodically review the procedures for operation and safety precautions.

1.1 System Overview

Indications for Use

The Digital Color Doppler Palm Ultrasound System is intended for diagnostic ultrasound imaging in B (2D), B/M, CFM, Combined (B+CFM), Pulsed Wave and Fusion Harmonic Imaging modes. It is indicated for Cardiac, Abdomen, Lung.

The Digital Color Doppler Palm Ultrasound System is intended for use in environments where healthcare is provided by healthcare professionals.

Contraindication

The system is NOT intended for ophthalmic use or any use that causes the acoustic beam to pass through the eye.

1.2 Contact Information

For additional information or assistance, please contact your local distributor or the appropriate support resource shown below:

CHISON website	www.chison.com
Service Support	CHISON Medical Technologies Co., Ltd. Tel: 0086-400-8878-020; 0086-0510-85311707 Fax: 0086-0510-85310726 E-mail: service@chison.com.cn
Placing an Order	CHISON Medical Technologies Co., Ltd. Tel: 0086-0510-8531-0593/0937 Fax: 0086-0510-85310726 Email: export@chison.com.cn
Manufacturer	CHISON Medical Technologies Co., Ltd. No.3 Changjiang South Road, Xinwu District, Wuxi, 214028 Jiangsu, P.R. China
US Agent	Mr. Marco Mu, 2219 Rimland Drive, Suite 301, Bellingham, Barkley Villiage Bellingham, Washington, 98226, UNITED STATES Phone: (702) 209-5185, Fax: 360-9253199 Email: us.agent@mid-link.net MID-LINK INTERNATIONAL CO., LTD



CAUTION

Federal law restricts the device to sale by or on the order of a licensed practitioner or therapist.

Chapter 2 System Safety

2.1 Safety Overview

This section discusses measures to ensure the safety of both the operator and patient. To ensure the safety of both operator and patient, please read the relevant details in this chapter carefully before operating this system. Disregarding the warnings or violation of relevant rules may result in personal injury or even loss of life for operator or patient.

Users should observe the following precautions:

- This system complies with Type BF general equipment, and the IEC standard. Please follow Chapter 1 “System Safety” in the user’s manual to use this system properly.
- Do not modify this system in any way. The system is prohibited to dismount. Necessary modifications must be made only by the manufacturer or its designated agents.
- This system has been fully adjusted at the factory. Do not adjust any fixed adjustable parts.
- In the event of a malfunction, turn off the system immediately and inform the manufacturer or its designated agents.
- Only connect this system, either electronically or mechanically, with devices that comply with the EN60601-1 standard. Recheck the leakage current and other safety performance indices of the entire system to avoid potential system damage caused by leakage from a current superposition.
- The system does not incorporate any specialized protective measures in the event it is configured with high-frequency operation devices. The operator should use caution in these types of applications.
- The system should be installed only by personnel authorized by the manufacturer. Do not attempt to install the system by yourself.
- Only an authorized service engineer may perform maintenance.
- Only a qualified operator, or someone under qualified supervision, should use the system.
- Do not use this system in the presence of flammable substances, otherwise an explosion may occur.
- Do not continuously scan the same part of a patient or expose the patient to prolonged scanning, otherwise it may harm the patient.
- When using the system for ultrasound testing, use only qualified ultrasound gel that complies with system standards.
- Do not unplug system when the system is in active operation. Always go to EXAM screen when you need to remove the system.
- To prevent from arm or neck injury, the operator should not stay at the same position for too long during patient scanning without taking break.
- Do not put the system near the liquid.



NOTE: To dispose of this product properly, Please call your local service department.

2.2 Electrical Safety

Type of protection against electric shock

Class I Equipment

CLASS I EQUIPMENT in which protection against electric shock does not rely on BASIC INSULATION only, but includes a protective earth ground. This additional safety precaution prevents exposed metal parts from becoming LIVE in the event of an insulation failure.



The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential ENVIRONMENT (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Degree of protection against electric shock

- **Type BF Applied part** (for Systems marked with BF symbol)

TYPE BF APPLIED PART providing a specified degree of protection against electric shock, with particular regard to allowable LEAKAGE CURRENT

BF: Isolation from ground; max. Patient leakage current: normal mode $\leq 100 \mu\text{A}$, single fault condition $\leq 500 \mu\text{A}$

Level of protection against harmful ingress of water

- The IP Classification of System is Ordinary Equipment, immersed part is IPX7, the other part is IPX1.

Safety level when used in the presence of FLAMMABLE ANAESTHETIC MIXED WITH AIR (or WITH OXYGEN or WITH NITROUS OXIDE):

The Equipment is not suitable for use in the environment with FLAMMABLE ANAESTHETIC MIXED WITH AIR (or WITH OXYGEN or WITH NITROUS OXIDE)

Conduction Interference



Image quality will effect by conducted disturbance. Please do not use the effected image. And put system in a simpler electromagnetic compatibility environment to start working, if conducted disturbance occurred, stop using it and change the position for acquiring a better image.

Mode of operation

- Continuous Operation

For maximum safety, always follow these guidelines:

- Do not remove the protective covers on the system. These covers protect users from hazardous voltages. A qualified electronic technician must make all internal replacements.
- Do not operate this system in the presence of flammable gases or anesthetics.

Notice upon Installation of Product

Separation distance and effect from fixed radio communications equipment: field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast transmitter cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ultrasound system is used exceeds the applicable RF compliance level as stated in the immunity declaration, the ultrasound system should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as re-orienting or relocating the ultrasound system or using an RF shielded examination room may be necessary.

- Use either power supply cords provided by or designated by CHISON. Never use any adaptor or converter to connect with a power source plug (e.g. three-prong- to-two-prong converter).
- Locate the equipment as far away as possible from other electronic equipment.
- Be sure to use only the cables provided by or designated by CHISON. Connect these cables following the installation procedures.



Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this medical system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Notice against User Modification










The user should never modify this product.

User modifications may cause degradation in Electrical Safety. Modification of the product includes changes in:

- Cables (length, material, wiring, etc.)
- System configuration/components

User modifications may cause degradation in EMC performance. Modification of the product includes changes in:


- Cables (length, material, wiring, etc.)
- System installation/layout

	This symbol is accompanied by the name and the address of the manufacturer and the manufacturing date of the device in the form YYYY-MM.
IPX7,IPX1	Protection against the effects of immersion
	This label indicates immersed part.
	Refer to instruction manual.
	This symbol indicates that in the united states of America, Federal law restricts the device to sale by or on the order of a licensed practitioner or therapist.
	WASTE OF ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE): This symbol is used for Environment Protection, it indicates that the waste of electrical and electronic equipment must not be disposed as unsorted waste and must be collected separately. Please contact your local Authority or distributor of the manufacturer for information concerning the decommissioning of your equipment.
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY: This symbol is accompanied by the name and the address of the authorized representative in the European Community.
	Scan the QR code to open the SonoEye website for more information.
	This symbol indicates the item is a medical device.
 (01) 0 6945121 40733 5 (11) 000000 (21) 123456789	This symbol indicates the UDI of the device, (01) is followed by the UDI-DI code of the device, (11) is followed by the manufacturing date of the device, (21) is followed by the serial number of the device.

2.3.2 Other Device Labels

The following table describes the purpose and location of safety labels and other important information provided on the equipment.

Table 2-1: Symbol Icons

Icon	Meaning
Identification and Rating Plate	<ul style="list-style-type: none"> • Manufacture's name • Serial number
	Type-BF applied part

2.4 Patient Environmental Devices

Front side:

- Power switch
- Freeze button, P button, + button, - button

Acceptable Devices

The Patient Environmental devices shown above are specified to be suitable for use within the PATIENT ENVIRONMENT.

Anyone using the equipment must be able to recognize the ESD symbol and understand how to take the necessary precautionary procedures, as described in the caution below:



- *DO NOT connect any device without approval by CHISON within the PATIENT ENVIRONMENT.*
- *DO NOT touch patient and devices without IEC/EN 60601-1 approval to avoid the leakage current risk within the PATIENT ENVIRONMENT.*

Unapproved Devices



- *DO NOT use unapproved devices.*
- *If devices are connected without the approval of CHISON, the warranty will be INVALID.*
- *The system can't be used with HF surgical equipment, otherwise the burns to patient may occur.*

Any device connected to this system must conform to one or more of the requirements listed below:

- *IEC standard or equivalent standards appropriate to devices.*



Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.



Unsafe operation or malfunction may result. Use only the accessories, options and supplies approved or recommended in these instructions for use.

Peripheral used in the patient environment

The system has been verified for overall safety, compatibility and compliance with the printer which is Mopria certified.



Printing quality may vary depending on the printer. If there is any printing quality problem, CHISON will not be responsible.

The system may also be used safely while connected to devices other than those recommended above if the devices and their specifications, installation, and interconnection with the system conform to the requirements of IEC/EN 60601-1-1.

The connection of equipment or transmission networks other than as specified in the user instructions can result in an electric shock hazard or equipment malfunction. Substitute or alternate equipment and connections require verification of compatibility and conformity to IEC/EN 60601-1-1 by the installer.

Equipment modifications and possible resulting malfunctions and electromagnetic interference are the responsibility of the owner.

 **WARNING**

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

2.5 Biological Safety

This product, as with all diagnostic ultrasound equipment, should be used only for valid reasons and should be used both for the shortest period of time and at the lowest power settings necessary (ALARA - As Low As Reasonably Achievable) to produce diagnostically acceptable images. The AIUM offers the following guidelines:

Clinical Safety Quoted from AIUM

Approved March 26, 1997

Diagnostic ultrasound has been in use since the late 1950s. Given its known benefits and recognized efficacy for medical diagnosis, including use during human pregnancy, the American Institute of Ultrasound in Medicine herein addresses the clinical safety of such use:

There are no confirmed biological effects on patients or instrument operators caused by exposures from present diagnostic ultrasound instruments. Although the possibility exists that such biological effects may be identified in the future, current data indicate that the benefits to patients of the prudent use of diagnostic ultrasound outweigh the risks, if any that may be present.

Heating: Elevating tissue temperature during obstetrical examinations creates medical concerns. At the embryo development stage, the rise in temperature and the length of time exposed to heat combine to determine potential detrimental effects. Exercise caution particularly during Doppler/Color exams. The Thermal Index (TI) provides a statistical estimate of the potential temperature elevation (in centigrade) of

tissue temperature. Three forms of TI are available: Soft Tissue Thermal Index (TIS), Bone Thermal Index (TIB).

Soft Tissue Thermal Index (TIS). Used when imaging soft tissue only, it provides an estimate of potential temperature increase in soft tissue.

Bone Thermal Index (TIB). Used when bone is near the focus of the image as in the third trimester OB examination, it provides an estimate of potential temperature increase in the bone or adjacent soft tissue.

Cavitation: Cavitation may occur when sound passes through an area that contains a cavity, such as a gas bubble or air pocket (in the lung or intestine, for example). During the process of cavitation, the sound wave may cause the bubble to contract or resonate. This oscillation may cause the bubbles to explode and damage the tissue. The Mechanical Index (MI) has been created to help users accurately evaluate the likelihood of cavitation and the related adverse effects.

MI recognizes the importance of non-thermal processes, cavitation in particular, and the Index is an attempt to indicate the probability that they might occur within the tissue.

2.6 Scanning Patients and Education

The Track-3 or IEC60601-2-37 output display standard allows users to share the responsibility for the safe use of this ultrasound system. Follow these usage guidelines for safe operation:

- In order to maintain proper cleanliness of the systems, always clean them between patients.
- Always use a disinfected sheath on all EV/ER systems during every exam.
- Continuously move the system, rather than staying in a single spot, to avoid elevated temperatures in one part of the patient's body.
- Move system away from the patient when not actively scanning.
- Understand the meaning of the TI, TIS, TIB and MI output display, as well as the relationship between these parameters and the thermal/cavitation bioeffect to the tissue.
- Expose the patient to only the very lowest practical transmit power levels for the shortest possible time to achieve a satisfactory diagnosis (ALARA - As Low As Reasonably Achievable).

2.6.1 Safe Scanning Guidelines

- Ultrasound should only be used for medical diagnosis and only by trained medical personnel.
- Diagnostic ultrasound procedures should be done only by personnel fully trained in the use of the equipment, in the interpretation of the results and images, and in the safe use of ultrasound (including education as to potential hazards).
- Operators should understand the likely influence of the machine controls, the operating mode (e.g. B-mode, color Doppler imaging or spectral Doppler) and system frequency on thermal and cavitation hazards.
- Select a low setting for each new patient. Output should only be increased during the examination if penetration is still required to achieve a satisfactory result, and after the Gain control has been moved to its maximum value.

- Maintain the shortest examination time necessary to produce a useful diagnostic result.
- Do not hold the system in a fixed position for any longer than is necessary. It should be removed from the patient whenever there is no need for real-time imaging or spectral Doppler acquisition. The frozen frame and Cine loop capabilities allow images to be reviewed and discussed without exposing the patient to continuous scanning.
- Take particular care to reduce output and minimize exposure time of an embryo or fetus when the temperature of the mother is already elevated.
- Take particular care to reduce the risk of thermal hazard during diagnostic ultrasound when exposing: an embryo less than eight weeks after gestation; or the head, brain or spine of any fetus or neonate.
- Operators should continually monitor the on-screen thermal index (TI) and mechanical index (MI) values and use control settings that keep these settings as low as possible while still achieving diagnostically useful results. In obstetric examinations, TIS (soft tissue thermal index) should be monitored during scans carried out in the first eight weeks after gestation, and TIB (bone thermal index) thereafter.

MI>0.3 There is a possibility of minor damage to neonatal lung or intestine. If such exposure is necessary, reduce the exposure time as much as possible.

MI>0.7 There is a risk of cavitation if an ultrasound contrast agent containing gas micro-spheres is being used. There is a theoretical risk of cavitation without the presence of ultrasound contrast agents. The risk increases with MI values above this threshold.

TI>0.7 The overall exposure time of an embryo or fetus should be restricted in accordance with **Table 2-2** below as a reference:

TI	Maximum exposure time (minutes)
0.7	60
1.0	30
1.5	15
2.0	4
2.5	1

Table 2-2 Maximum recommended exposure times for an embryo or fetus

- Non-diagnostic use of ultrasound equipment is not generally recommended. Examples of non-diagnostic uses of ultrasound equipment include repeated scans for operator training, equipment demonstration using normal subjects, and the production of souvenir pictures or videos of a fetus. For equipment of which the safety indices are displayed over their full range of values, the TI should always be less than 0.5 and the MI should always be less than 0.3. Avoid frequent repeated exposure of any subject. Scans in the first trimester of pregnancy should not be carried out for the sole purpose of producing souvenir videos or photographs, nor should their production involve increasing the exposure levels or extending the scan times beyond those needed for clinical purposes.
- Diagnostic ultrasound has the potential for both false positive and false negative results. Misdiagnosis is far more dangerous than any effect that might result from the ultrasound exposure. Therefore,

diagnostic ultrasound system should be performed only by those with sufficient training and education.

2.6.2 Understanding the MI/TI Display

Track-3 follows the Output Display Standard for systems that include fetal Doppler applications. The acoustic output will not be evaluated on an application-specific basis, but the **global maximum de-rated Ispta** must be $\leq 720 \text{ mW/cm}^2$ and either the **global maximum MI** must be ≤ 1.9 or the **global maximum de-rated Ispta** must be $\leq 190 \text{ W/cm}^2$. An exception is for ophthalmic use, in which case the $\text{TI} = \max(\text{TIS}_{\text{as}})$ is not to exceed 1.0; $\text{Ispta.3} \leq 50 \text{ mW/cm}^2$, and $\text{MI} \leq 0.23$. **Track-3** gives the user the freedom to increase the output acoustic power for a specific exam, and still limit output acoustic power within the **global maximum de-rated Ispta** $\leq 720 \text{ mW/cm}^2$ under an Output Display Standard.

For any diagnostic ultrasonic systems, Track-3 provides an Output Indices Display Standard. The diagnostic ultrasound systems and its operator's manual contain the information regarding an ALARA (As Low As Reasonably Achievable) education program for the clinical end-user and the acoustic output indices, MI and TI. The MI describes the likelihood of cavitation, and the TI offers the predicted maximum temperature rise in tissue as a result of the diagnostic examination. In general, a temperature increase of 2.5°C must be present consistently at one spot for 2 hours to cause fetal abnormalities. Avoiding a local temperature rise above 1°C should ensure that no thermally induced biologic effect occurs. When referring to the TI for potential thermal effect, a TI equal to 1 does not mean the temperature will rise 1 degree C. It only means an increased potential for thermal effects can be expected as the TI increases. A high index does not mean that bioeffects are occurring, but only that the potential exists and there is no consideration in the TI for the scan duration, so minimizing the overall scan time will reduce the potential for effects. These operator control and display features shift the safety responsibility from the manufacturer to the user. So it is very important to have the Ultrasound systems display the acoustic output indices correctly and the education of the user to interpret the value appropriately.

RF: (De-rating factor)

In Situ intensity and pressure cannot currently be measured. Therefore, the acoustic power measurement is normally done in the water tank, and when soft tissue replaces water along the ultrasound path, a decrease in intensity is expected. The fractional reduction in intensity caused by attenuation is denoted by the de-rating factor (RF),

$$\text{RF} = 10^{(-0.1 a f z)}$$

Where a is the attenuation coefficient in $\text{dB cm}^{-1} \text{ MHz}^{-1}$, f is the transducer center frequency, and z is the distance along the beam axis between the source and the point of interest.

De-rating factor RF for the various distances and frequencies with attenuation coefficient 0.3 dB

cm-1 MHz-1 in homogeneous soft tissue is listed in the following table. An example is if the user uses 7.5MHz frequency, the power will be attenuated by .0750 at 5cm, or $0.3 \times 7.5 \times 5 = -11.25\text{dB}$. The De-rated Intensity is also referred to as '.3' at the end (e.g. Ispta.3).

	Distance		Frequency (MHz)		
	(cm)	1	3	5	7.5
1	0.9332	0.8128	0.7080	0.5957	
2	0.8710	0.6607	0.5012	0.3548	
3	0.8128	0.5370	0.3548	0.2113	
4	0.7586	0.4365	0.2512	0.1259	
5	0.7080	0.3548	0.1778	0.0750	
6	0.6607	0.2884	0.1259	0.0447	
7	0.6166	0.2344	0.0891	0.0266	
8	0.5754	0.1903	0.0631	0.0158	

$I' = I \cdot R^2$ Where I' is the intensity in soft tissue, I is the time-averaged intensity measured in water.

Tissue Model:

Tissue temperature elevation depends on power, tissue type, beam width, and scanning mode. Six models Tissue temperature elevation depends on power, tissue type, beam width, and scanning mode. Six models are developed to mimic possible clinical situations.

	Thermal Models	Composition	Mode	Specification	Application
1	TIS	Soft tissue	Unscanned	Large aperture (>1cm ²)	Liver PW
2	TIS	Soft tissue	Unscanned	Small aperture (<1cm ²)	Pencil System
3	TIS	Soft tissue	Scanned	Evaluated at surface	Breast color
4	TIB	Soft tissue and bone	Scanned	Soft tissue at surface	Muscle color
5	TIB	Soft tissue and bone	Unscanned	Bone at focus	Fetus head PW

Soft tissue:

Describes low fat content tissue that does not contain calcifications or large gas-filled spaces.

Scanned: (auto-scan)

Refers to the steering of successive burst through the field of view, e.g. B and CFM mode.

Unscanned:

Emission of ultrasonic pulses occurs along a single line of sight and is unchanged until the transducer is moved to a new position. For instance, the PW mode.

TI:

TI is defined as the ratio of the In Situ acoustic power ($W.3$) to the acoustic power required to raise tissue temperature by 1°C (W_{deg}), $TI = W.3/W_{deg}$.

Three TIs corresponding to soft tissue (TIS) for abdominal; bone (TIB) for fetal and neonatal cephalic; have been developed for applications in different exams.

An estimate of the acoustic power in milliwatts necessary to produce a 1°C temperature elevation in soft tissue is:

$$W_{deg} = 210/f_c, \quad \text{for model 1 to 4, where } f_c \text{ is the center frequency in MHz.}$$

$$W_{deg} = 40 K D^2 \quad \text{for model 5 and 6, where } K \text{ (beamshape factor) is } 1.0, D \text{ is the aperture diameter in cm at the depth of interest.}$$

MI:

Cavitation is more likely to occur at high pressures and low frequencies in pulse ultrasound wave in the tissue, which contains the bubble or air pocket (for instance, the lung, intestine, or scan with gas contrast agents). The threshold under optimum conditions of pulsed ultrasound is predicted by the ration of the peak pressure to the square root of the frequency.

$$MI = Pr' / \text{sqrt}(f_c)$$

Pr' is the de-rated (0.3) peak rare-fractional pressure in Mpa at the point where PII is the maximum, and f_c is the center frequency in MHz. PII is the Pulse Intensity Integral that the total energy per unit area carried by the wave during the time duration of the pulse. The peak rare- fractional pressure is measured in hydrophone maximum negative voltage normalized by the hydrophone calibration parameter.

Display Guideline:

For different operation modes, different indices must be displayed. However, only one index needs to be shown at a time. Display is not required if maximum MI is less than 1.0 for any setting of the operating mode, or if maximum TI is less than 1.0 for any setting of the operating mode. For TI, if the TIS and TIB are both greater than 1.0, the scanners need not be capable of displaying both indices simultaneously. If the index falls below 0.4, no display is needed.

Display and Report in Different Mode

Located on the upper middle section of the system display monitor, the acoustic output display provides the operator with real-time indication of acoustic levels being generated by the system.

Only display and report TIS or TIB and start from 0.4 if maximum TI > 1.0, display in increments of 0.2 for values of indices of 2.0 or less, and 0.5 for values of indices greater than 2.0.

Below is a simple guideline for the user when TI exceeds one limit exposure time to $4(6-TI)$ minutes based

on the 'National Council on Radiation Protection. Exposure Criteria for Medical Diagnostic Ultrasound: I. Criteria Based on Thermal Mechanisms. Report No.113 1992'.

Operator Control Features:

The user should be aware that certain operator controls may affect the acoustic output. It is recommended to use the default (or lowest) output power setting and compensate using Gain control to acquire an image. Other than the output power setting in the soft-menu, which has the most direct impact on the power; the PRF, image sector size, frame rate, depth, and focal position also slightly affect the output power. The default setting is normally around 70% of the allowable power depending on the exam application mode.

Controls Affecting Acoustic Output

The potential for producing mechanical bioeffects (MI) or thermal bioeffects (TI) can be influenced by certain controls.

Direct: The Acoustic Output control has the most significant effect on Acoustic Output.

Indirect: Indirect effects may occur when adjusting controls. Controls that can influence MI and TI are detailed under the Bioeffects portion of each control in the Optimizing the Image chapter.

Always observe the Acoustic Output display for possible effects.

Best practices while scanning

HINTS: Raise the Acoustic Output only after attempting image optimization with controls that have no effect on Acoustic Output, such as Gain and STC.



WARNING: *Be sure to have read and understood control explanations for each mode used before attempting to adjust the Acoustic Output control or any control that can effect Acoustic Output.*

Use the minimum necessary acoustic output to get the best diagnostic image or measurement during an examination. Begin the exam with the system that provides an optimum focal depth and penetration.

Acoustic Output Default Levels

In order to assure that an exam does not start at a high output level, the system initiates scanning at a reduced default output level. This reduced level is preset programmable and depends upon the exam icon and system selected. It takes effect when the system is powered on or New Patient is selected. To modify acoustic output, adjust the Power Output level on the Soft Menu.

2.7 Device instructions of cybersecurity controls

2.7.1 Software Integrity Control

The following controls are in place to assure that the device software will maintain its integrity from the point of origin to the point at which that device leaves the control of the manufacturer.

- The system goes through a security self-test at startup, some key files for security checks to ensure that the system starts normally.
- The system supports the recovery function. This function will restore the system to the original state. The following controls are in place to assure that the device software will be protected from malware from the point of origin to the point at which that device leaves the manufacturer's control point.
- Installing OS Updates Immediately
Whenever an update appears on the screen, consider installing an operating system update immediately to avoid virus and malware attacks.
- The device software doesn't provide any entrance for executing third party application.
- Never trust an unknown computer
The device connect computer with USB port. Do not plug the device into an unknown computer.
- Install anti-virus software and use firewall and scan device regularly with antivirus software.
- Keep network safe
Set a strong network access password, and do not open public Wi-Fi connections, use WPA, WPA2 encryption or the latest WPA3 encryption. Don't open email attachments from unknown people or companies, don't click links in unsolicited emails, do not download suspicious apps, etc.

2.7.2 Device instructions of cybersecurity controls

The following are device instructions for use related to recommended cybersecurity controls appropriate for the intended use environment.

- Install anti-virus software and use firewall before connecting the device with wifi.
- Install anti-virus software and use firewall before connecting the device with DICOM server.
- Limit access to device software through the authentication of users by user name and password. Ownership of a device is assigned to one user at a time.
- To protect the patient information, the system should hide the critical patient information when exporting the image and cine. And the hidden function is configurable.
- Image and cine data contains no patient or user-identifying information. If want to encrypt this data, please connect to a network that uses an encryption protocol.
- The software uses private custom format to store patient data and cannot be recognized by general tools.
- Connect the network only to run the DICOM functions or transmission functions. Otherwise disconnect the network.
- When connecting the device with wifi, use a network that supports Wi-Fi 802.11n. We recommend that secure this network using WPA (Wi-Fi Protected Access) or WPA2 (Wi-Fi Protected Access II) as security protocol.

2.8 Instructions to the intruded device

The following are instructions to protect your device prior to the device leaving your control.

- Back up the patient data.

- Verify the patient data on the device with anti-virus software.
- Recovery the system to the original state.
- Recovery the patient data to the system.

CHISON

Value Beyond Imaging



CHISON Medical Technologies Co., Ltd.

No.3 Changjiang South Road, Xinwu District, Wuxi, 214028 Jiangsu, P.R. China