



FM-200

Fetal Doppler



FM-200 Fetal Doppler

Water-proof Probe



Back light LCD



Operation buttons



Carry bag



- Features:**
- Handheld & elegant design;
 - Back light LCD;
 - Display of FHR, storage, work mode, battery status, signal indication and probe frequency;
 - Three work modes: real time, average and manual;
 - FHR alarm limits is adjustable and can be saved;
 - Low power consumption design;
 - Unique heart rate algorithm enable FHR more accurate and reliable;
 - Buckle design for easy take;
 - Automatically power off while no signal within 2 minutes;
 - Inter-changeable water-proof probe of 2MHz or 3MHz.

Options: 3MHz probe, Carry bag, Earphone, rechargeable battery.

Specifications:

Ultrasonic Work Frequency: 2MHz, 3MHz(Optional);	Power Requirement: 2 x 1.5V (AA) Alkaline Battery or Ni-MH Rechargeable Battery;
Heart Rate Range: 50bpm~240bpm;	Dimension: 130mm×89mm×34mm;
Alarm Limit: Adjustable lower limit range: 50bpm~120 bpm;	Net Weight: 450g(including battery).
Adjustable upper limit range:160bpm~240 bpm;	

Rabat : 2, rue Tissa, Av Allal ben abdellah, 10 000 Rabat Centre. **Tél :** 05 37 70 06 04 / 05 37 70 48 27 - **Fax :** 05 37 70 51 06
Casablanca : 52, Boulevard Abdelmoumen Résidence Al Manar. **Tél :** 05 22 25 41 03 / 05 22 25 09 10 - **Fax :** 05 22 99 35 89
Fes : 78, Quartier EL WAFA, Résidence Safae. **Tél :** 05 35 61 65 53 / 54 - **Fax :** 05 35 61 65 52
Marrakech : 32 Av. Ibn Sina Lotissement Boumessmar. **Tél :** 05 24 30 71 94 / 05 24 30 74 04 - **Fax :** 05 24 30 71 81

Email : topmedical@menara.ma

Site web : www.topmedical.ma

Facebook : Top Médical Maroc - La Page Officielle.