

Fingertip Pulse Oximeter



Top Médical®
Equipements de santé et de bien-être depuis 1979

Fingertip Pulse Oximeter



BP-10M

Features:

- OLED two-color display;
- Accurately measures SpO2 value and Pulse Rate value;
- Innovative 2 directions display;
- Automatic start measuring after putting finger into tuber cushion;
- Power off automatically when no signal;
- Audio & visible alarm function;
- Low battery voltage indication;
- Low power consumption, 30 hours battery life;
- 1/3 of business card size.




BP-10P

Features:

- OLED color display;
- Accurately measures SpO2 value and Pulse Rate value;
- 4-direction display automatically;
- Automatic start measuring after putting finger into tuber cushion;
- Power off automatically when no signal;
- Audio & visible alarm function;
- Low battery voltage indication;
- Low power consumption, 30 hours battery life;
- 1/3 of business card size.

Specifications:

Picture			
Model		BP-10M	BP-10P
Display		OLED	OLED
Display Direction		2	4
Waveform		YES	
Audio Alarm		YES	
SPO2	Measurement Range	35%-99%	
	Accuracy	±3%(70%-75%), ±2%(75%-99%)	
	Resolution	1%	
PR	Measurement Range	30-240BPM	
	Accuracy	±2BPM	
	Resolution	1BPM	
PI	Measurement Range	0.2%-20%	
	Accuracy	±0.1%(during 0.2%-2%), ±1%(during 2%-10%)	
	Resolution	0.1%	
Low Perfusion		≤0.6%	
Alarm Indication	SPO2		< 90%
	PR	Upper Limit	120BPM
		Lower Limit	50BPM
Power Supply Management		Power off automatically (no signal for more than 8 seconds)	
Power Requirements		AAA size battery(1.5V)×2	
Power Consumption		<40mA	
Battery Life		30 hrs	
Dimension		66×36×33mm(L×W×H)	
Net Weight Including Battery		60g	

Rabat : 2, rue Tissa, Av Allal ben abdellah, 10 000 Rabat Centre. **Tél :** 05 37 70 06 04 / 05 37 70 48 27 - **Fax :** 05 37 70 51 06

Casablanca : 52, Boulevard Abdelmoumen Résidence Al Manar. **Tél :** 05 22 25 41 03 / 05 22 25 09 10 - **Fax :** 05 22 99 35 89

Fes : 78, Quartier EL WAFA, Résidence Safae. **Tél :** 05 35 61 65 53 / 54 - **Fax :** 05 35 61 65 52

Marrakech : 32 Av. Ibn Sina Lotissement Boumessmar. **Tél :** 05 24 30 71 94 / 05 24 30 74 04 - **Fax :** 05 24 30 71 81

Email : topmedical@menara.ma

Site web : www.topmedical.ma

Facebook : Top Médical Maroc - La Page Officielle.