

Top Médical

Equipements de santé et de bien être depuis 1979

iC 60

Fetal & Maternal
Monitor



iC 60 Fetal & Maternal Monitor



- Features:**
- Light & compact, user friendly design;
 - 7 inches high resolution foldable TFT color screen;
 - Basic parameters: FHR, TOCO, AFM, M HR, NIBP, SPO2;
 - Accuracy & reliability DSP technology for real time measurement and identification;
 - Patient event marker and clinical event marking button mark clinical events separately;
 - Provide convenience and reliability during antepartum, intrapartum and postpartum monitoring;
 - 12 hours CTG storage, fast review and print, optional to upgrade for 24 hours;
 - Unique design for cable of transducers to avoid intertwisting;
 - Alarm system: time-over alarm, non-paper alarm, out of range FHR alarm;
 - Automatically adjustable many kinds of 110mm thermal paper with different specifications;
 - Wire or wireless connect to central monitoring system;
 - Low acoustic power, safer to the fetus.

- Options:**
- Twins FHR measurement;
 - Build-in network/PC connecting interface;
 - FAS(fetal acoustic stimulator);
 - Wireless network.

Specification:

FHR

Ultrasound frequency: 1.0 MHz
 Ultrasound strength: <5mW/cm²
 FHR Range: 50-210 bpm
 Accuracy: ±2 bpm

TOCO

TOCO Range: 0-100 units
 Non-linear tolerance: <±10%

Auto Fetal Movement

Range: 0-100

Maternal SpO2

Measurement Range: 70% ~ 100%
 Accuracy: ±3%

HR

Measurement Range: 30-240 bpm
 Accuracy: ±2 bpm

Maternal NIBP

Systolic pressure: 50-240mmHg
 Mean pressure: 15-180mmHg
 Diastolic pressure: 25-200mmHg

Dimension: 295mm×240mm×73mm(L×W×H)
 Net Weight: 1.75 Kg

Rabat : 2, rue Tissa, Av Allal ben abdellah, 10 000 Rabat Centre. **Tél :** 05 37 70 06 04 / 05 37 70 48 27 - **Fax :** 05 37 70 51 06
Casablanca : 52, Boulevard Abdelmoumen Résidence Al Manar. **Tél :** 05 22 25 41 03 / 05 22 25 09 10 - **Fax :** 05 22 99 35 89
Fes : 78, Quartier EL Wafa, Résidence Safae. **Tél :** 05 35 61 65 53 / 54 - **Fax :** 05 35 61 65 52
Marrakech : 32 Av. Ibn Sina Lotissement Boumessmar. **Tél :** 05 24 30 71 94 / 05 24 30 74 04 - **Fax :** 05 24 30 71 81

Email : topmedical@menara.ma

Site web : www.topmedical.ma

Facebook : Top Médical Maroc - La Page Officielle.