



## NEW POCKET PHYSIO PRO

Portable professional electrostimulator for physiotherapy; software with emission of all electrotherapy currents for analgesic and anti-inflammatory applications, post-injury rehabilitation of healthy and denervated muscles, urogynecological therapy. It includes pre-set programmes for the most common applications and it is possible to set new parameters and memorize new programmes. The 2 independent channels allow a differentiated stimulation in the different muscles and the interactive display allows the contemporary selection of 2 different programmes (for example in the case of the simultaneous treatment of two patients who have different needs). The software includes 135 pre-set programmes and 30 free. Double feeding: rechargeable battery and power mains.



### Programmes / Applications

Anti-inflammatory and Decontracting Effect	Post-injury Functional Rehabilitation
Trophic Effect for Traumas with Edemas	Relaxing Effect for Tired Muscles
Capillarisation	Stimulation of Denervated Muscles
Cervical Pain	Stimulation of Normoinnervated Muscles
Sprain	Muscle Strain
Acute and Chronic Pains	Strength Development
Epicondylitis	Tendonitis
Knee Pain	Toning
Ionophoresis for Pains and Inflammations	Analgesic Urogynecological Applications
Lumbago	Hyperhidrosis Applications
Muscular Warm-up	Carpal Tunnel Syndrome
Recovery from Muscular Fatigue	Vasodilation and Hyperemising Effect

### Technical Specifications

<b>Constant Current Generator</b>	yes, 1	<b>Standard Kit</b>
<b>Power Supply</b>	rechargeable battery / power mains	Bag
<b>Wave Forms</b>	tens, diadynamic, bipolar interferential, rectangular biphasic, faradic, kotz, ionophoresis, rectangular monophasic, triangular, exponential, uro	1 set of self-adhesive plug electrodes (4 50x50mm, 2 50x100mm)
<b>No. of Output Channels</b>	2 independent	4 conductive rubber electrodes (2 50x50mm, 2 60x80mm)
<b>No. of Programmes</b>	135 pre-set + 30 free	4 electrode sponges (2 50x50mm, 2 60x80mm)
<b>Current Intensity (mA)</b>	0-120	2 elastic bands for electrodes
<b>Frequency (Hz)</b>	1-160	Rechargeable battery
<b>Emission</b>	continuous/modulated	Power supply
<b>Impulse duration</b>	50 µs - 1 sec	Output cables
<b>Dimensions and weight</b>	18*11*4cm, 400 gr only device, 1,5 kg with accessories	2 double-end cables
		User's manual with photos of the applications