



NEW POCKET PHYSIO IONOTENS

The Pocket Physio Line includes a new electrostimulator specific for analgesic and anti-inflammatory applications: the IONOTENS model. Its software includes 33 pre-set programmes for the most common applications and 4 free programmes for setting personalised programmes. The double feeding: rechargeable battery and power mains. It is so easy to use that the patient can make applications alone at home. It issues analgesic tens, ionophoresis, iontoforesis and anti-inflammatory currents



Programmes / Applications

Acute and Chronic Pains	Strains
Sprains	Tears
Knee Pain	Contracture
Cervical Pain	Bruise
Lumbago	Ionophoresis for pains
Epicondylitis	Ionophoresis for acute inflammations
Tendonitis	Ionophoresis for chronic inflammations
Carpal Tunnel Syndrome	

Technical Specifications

Constant Current Generator	yes, 1
Power Supply	rechargeable battery / power mains
Wave Forms	tens, ionophoresis
No. of Output Channels	2 independent
No. of Programmes	34 pre-set + 4 free
Current Intensity (mA)	0-120
Frequency (Hz)	1-160
Emission	continuous/modulated
Impulse duration	50 μ s - 1 sec
Dimensions and weight	18*11*4cm, 400 gr only device, 1,5 kg with accessories

Standard Kit

Bag
1 set of self-adhesive plug electrodes (4 50x50mm, 2 50x100mm)
4 conductive rubber electrodes (2 50x50mm, 2 60x80mm)
4 electrode sponges (2 50x50mm, 2 60x80mm)
2 elastic bands for electrodes
Rechargeable battery
Power supply
Output cables
2 double-end cables
User's manual with photos of the applications